MARK · VU

It is easy to observe a skin problems and compare before & after due to rotate plate makes to easy taking a front, left, and right side of facial skin








Left 30°

MERIT








Capture button Capture button for getting a great reproducibility



Power & USB ·Connected to PC by USB ·Power : DC 19V 3.42A

LIFETIME

It is possible to take 15 million shots through LED light without replaced light sources.



End-of-life for LEDs

Regulatory agencies are specifying

ENERGY STAR SSL Criteria V.1.1 specifies minimun lifetime Lumen maintenance, L70:25,000 hours for residential; 35,000 hours for commercial

·Mark-Vu LED lifetime = More than 150,000 shots available ·Other companies Flash lifetime = 3000~5000 shots available *When shooting 150,000 shots, cost of consuming more than 20,000 dollar.





MARK · Vu

By applying Canon DSLR camera and lens, it is possible to make a high quality image that can be used in labs and dermatology.

Certification & Patent

10-2010487 Certificate of Patent No. Design Registration No. 30-0935733 40-1334710 Trademark Registration No.



GHOST MIRROR

Mirror type facial skin analyzer that can directly observe the skin problems through the mirror by yourself.

Certification & Patent

Certificate of Patent No. 10-1981997 30-0920516 Design Registration No. Trademark Registration No. 40-1296863

LICENSE









PSPLIS COLID
S. Monorow, Younghong o, O;
Sopulitie of Korea
Sopulitie of Korea
Sopulities (States of Manufacture)

















Mark-Vu

在 图 人 SECSESCO



Office

Factory

102-204, 88, Sinwon-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, Republic of Korea, 16681 Tel. +82-(0)31-695-6916 Fax. +82-(0)31-695-6919 Email. psiplus@psiplus.co.kr

102-408, 88, Sinwon-ro, Yeongtong-gu, Suwon-si, Gyeonggi-do, Republic of Korea, 16681 Tel. +82-(0)31-695-6917 Fax. +82-(0)31-695-6919 Web. www.psiplus.co.kr



Professional Facial Skin Analyzer





MARK · VU

Professional Facial Skin Analyzer

Mark-Vu, a professional skin analyzer that uses all light sources as LED, delivers the best image.

The filter on / off system using a motor makes a clear image from ultraviolet light and normal light. By using curtains and a rotating plate, it is possible to make a clearer image. It is also easy to use by using the capture button at the top of the equipment.



LIGHT SOURCE

Patent No. of light sources 10-2010487 Patent No. of Hidden Wrinkles 10-1981997

Since all four different kinds of lights sources are applied with LED, accurate before & after images can be seen.





Pore Wrinkles Skin tone Dark-Circle




Hidden Wrinkles Skin Texture Radiance



Polarized

Hidden Wrinkles Fine Wrinkles Skin Texture



Melanin Sebum Porphyin

HIGH QUALITY

Easy to observe skin problems clearly, as the images are being captured by light resolution DSLR camera




DSLR Camera capture images

Detail images

FUNCTION

There are added comport function as like Sorting as date, same light sources, Search a client and etc.







Button for Search a client :

You don't need to move to STEP 1 for search a client due to there is button for search a client at STEP 3.






Button for making same light source Every light sources will be changed to same light source as you choose.









RESULT

It is easy to consult with the 12 kinds of skin problems identified and customized description according to the customer's skin condition



















COMPARISON

By comparing before & after images, improvement of skin issues can be easily acknowledged.



