

### **Clinical Picture** Aesthetic

Wrinkle



Before & After

# **Technical Specification**

Standa	Standard type	
Freq	uency	4MHz,7MHz
Max E	Energy	3.0J/cm2
Cart	ridge	7MHz:1.5mm, 4MHz:10mm,1
Spa	ncing	1.0~2.0mm(0.1
Powe	r Level	Adjustable 1~1
Ler	ngth	Full Length:25
Spot	t Size	<3mm
	CD	10.2Inchs True
Cart Spa Powe Ler Spot	ridge acing r Level ngth t Size	7MHz:1.5mm, 4MHz:10mm, 1.0~2.0mm(0 Adjustable 1~ Full Length:25 <3mm

### Indications

• Smart mode(Preset mode) at each treatment

Skin iftina

- Real Penetration of Depth
- Skin tightening &Lifting
- Short treatment time
- Easy to use small area by pen type handpiece
- Contraction of dermis & SMAS layer by HIFU
- Non Surgical face lifting and tightening
- Face&Body(fat) treatment by one device

### **Technical Specification**

Pen type	
Frequency	7MHz
Max Energy	3.0J/cm2
Cartridge	7MHz:1.5mm, 3
Power Level	Adjustable 1~1
Speed	1-4shot
Spot Size	<3mm

### Functions

3.0mm, 4.5mm-For Face I 3mm(option)-For Body & Fat

1mm step)-For Fat

OLevel

Color Touch

Functions

3.0mm

10Level

ULTRA POWERFUL HIFU SYSTEM **BU:TERA** 





# **HIFU** BU:TERA

From the deep SMAS of skin to the dermis and the epidermis

# **VARIOUS TREATMENTS** AT ONCE!!



Use freely from small to large area.

# Cartridge Type

Cartridge type	D1.5	D3.0
Frequency	7MHz	7MHz
Depth of the skin	1.5mm	3.0mm
Pen type	1.5mm	3.0mm
	Standard	Standard

### Pen type cartridges For Face:1.5mm,3.0mm



P3.0

Standard cartridges For Face:1.0mm.3.0mm,4.5mm For Body:10mm,13mm



D4.5	F10	F15		
7MHz	4MHz	4MHz		
4.5mm	10mm (Fat)	15mm (Fat)		
Standard	Optional	Optional		



--> Standard Handpiece



# Smart mode installed for face and body







