

IDS SERIES



IDS-1000



IDS-1000 PLUS



Your Dermoscopy become now compact and Light

- Small enough for your pocket or hand bag
- Light enough for hanging on your neck or carry at hand the whole day
- Perfect for laypeople who want observe their skin changes personally
- Reasonably priced

IDS-1000 series are reasonably priced, but the functions are still the same

- 10 x magnification
- Contact & non-contact dermoscopy with polarized light
- Contact & non-contact dermoscopy with non-polarized light
- Smartphone & camera compatibility (IDS-1000 Plus only)

Excellent Optical and Illumination Quality

Dermoscopy requires

- ✓ a high-quality magnifying lens and likewise
- ✓ high quality light sources

that can help with examining skin structures and patterns.

As the IDS-1100, IDS-1000 has also developed by optical experts and we guarantee the excellent optical & light performance of IDS-1000 Series.

IDS-1000 Series are lightly portable and wearable Dermatoscopes made for anytime and anywhere use.

Compared to IDS-1100, IDS-1000 Series are small, light and installed slightly smaller lens, but the optical and illumination quality are the same, even the depth of field is deeper than IDS-1100.

	IDS-1000 Series	IDS-1100 Series
Polarization	Cross & Parallel	
Field of view	18 mm	23 mm
Magnification	10 x	
Depth of field	1.5 mm	1.3 mm
Resolution	95 LP	
Contact & non-contact examination	Yes	



IDS-1000 Series models are reasonably priced but the main quality and functions are still the same.

Polarized and non-polarized mode change

The Polarized and non-polarized mode are manually to switch by turning the cone hood.

- Polarization mode: 3 o'clock
- Non-polarization mode: 12 o'clock



Silicon Sleeve with lanyard

Your Dermatoscope will never be lost with the sleeve and lanyard.



Cone hood with contact plate

Differently from IDS-1100 Dermatoscope, the contact plates of IDS-1000 models are bonded to the cone hood and there are three different types of cone hood available.

- Cone hood (open top : with no contact plate) : for non-contact examination
- Cone hood contact plate (with no scale)
- Cone hood scale contact plate (with scale)



Cone hood (open top)



Cone hood contact plate (no scale)



Cone hood scale contact plate

Smartphone and camera compatibility

The only difference between IDS-1000 and IDS-1000Plus is the smartphone and camera compatibility.

IDS-1000Plus is compatible with smartphones, mirrorless cameras as well as compact cameras.

Attach your smartphone to the IDS-1000Plus, observe and take images of common skin lesions. If you need more detailed observation, use a camera to obtain dedicated and precise images.




Mirrorless camera adapter
for Canon, Sony, Nikon



Universal smartphone adapter



Smartphone adapter



Adhesive adapter unit for
smartphone & tablet pc

IDS-1000



Technical Specification

Polarization	Cross & Parallel
Aperture	18 mm
Magnification	10 x
Depth of field	1.5
Resolution	95 LP
LED	8 White LEDs
Focus Mode	Fixed
Battery	Disposable Lithium 6V 2CR5
Continuous duty time	8 hours
Camera & Smartphone Compatibility	X
Contact & Non-Contact Examination	Yes
Net Dimension / Weight	51 x 112 x 23 (mm) / 100g
Package Dimension / Weight	70 x 170 x 60 (mm) / 330g
Basic Components	<ul style="list-style-type: none">① IDS-1000② Cone hood (open top)③ Cone hood contact plate (no scale)④ Silicon sleeve with lanyard

IDS-1000 PLUS



Technical Specification

Polarization	Cross & Parallel
Aperture	18 mm
Magnification	10 x
Depth of field	1.5
Resolution	95 LP
LED	8 White LEDs
Focus Mode	Fixed
Battery	Disposable Lithium 6V 2CR5
Continuous duty time	8 hours
Camera & Smartphone Compatibility	Yes
Contact & Non-Contact Examination	Yes
Net Dimension / Weight	51 x 112 x 23 (mm) / 110g
Package Dimension / Weight	70 x 170 x 60 (mm) / 330g
Basic Components	<ul style="list-style-type: none">① IDS-1000Plus② Cone hood (open top)③ Cone hood scale contact plate④ Silicon sleeve with lanyard

ILLUCO 

102-304 SK Ventium, #166 Gosan-ro, Gunpo-si, 15850

Gyeonggi-do, South Korea

Tel : +82 31 429 8825

Int'l : +82 70 4228 2351

Fax : +82 31 429 8826

E-mail : info@illuco.co.kr

Web : www.illuco.co.kr