

## Digital Infusion Pump

# DermaShine PRO

CE MARK

Medical Device

### PURPOSE

- To inject HA into the skin surface using the multi needles
  - **DermaShine PRO** : Automated mesotherapy injector
  - **9 Pin Needle** : Sterile needle to inject drugs into the skin (disposable)

### EFFECTS

- Hydrated skin by injecting HA
- Radiant and taut skin, Pore reduction, Wrinkle improvement
- Less pain with shorter treatment time

### SPECIFICATIONS

Model	Panace-DS-51
Size (WxLxH)	379.7x328.0x178.4mm
Weight	3.4kg
Power rate	AC100-240V / 50/60Hz
Modes	Auto-Sensing, Auto-Dose, Dose, Continuous
Needle	Multi Needle (9 Pin)
Components	Injector, Suction Tube, Power Code, Adapter, Foot Switch, Stand(Optional)



### 9 Pin Multi Needle



- 9 Pin Multi Needle (32G / 34G) minimizes the patient's pain and enables them to return to their daily lives after the treatment.
- Able to customize the injection depth(0~2mm) according to the skin condition of the patients.